


Hot Air Oven

Brand Name: Advance Science-AS



CONSTRUCTION

- Double walled inside chamber made of Aluminium or Stainless Steel 304 qty.
- Outer body made of Mild Steel and painted in Epoxy Powder Coating.
- Insulation by Mineral wool 75mm.
- The door has synthetic neoprene rubber gasket.
- Supplied with 2/3 removable wire-mess shelves.
- Supplied with / without air circulating fan.

TEMPERATURE

- Thermostatic control the temperature range of 50°C to 250°C ±2°C.
- With coils heaters placed at the bottom
- With three controls switch (Rotary Switch) low, medium & high for heater/fan motor speed.
- Graduated Temperature control knob.
- Ventilators on both sides.
- Work on 220/230V AC, 50 Hz but without thermometer

SIZE INSIDE CHAMBER (ALUMINUM)

MODEL	SIZE (MM) W H D	LTRS	Price Without Fan	Price with Air Circulating Fan
HAO-1	300 x 300 x 300	28.0	N/A	N/A
HAO -2	355 x 355 x 355	45.0	N/A	N/A
HAO -3	455 x 455 x 455	95.0	N/A	N/A
HAO -4	455 x 455 x 605	125.0	N/A	N/A
HAO -5	605 x 605 x 605	224.0	N/A	N/A
HAO -6	605 x 455 x 910	252.0	N/A	N/A
HAO -7	605 x 605 x 910	336.0	N/A	N/A

SIZE INSIDE CHAMBER (STAINLESS STEEL)

MODEL	SIZE (MM) W H D	LTRS	Price Without Fan	Price with Air Circulating Fan
HAO -1	300 x 300 x 300	28.0	N/A	N/A
HAO -2	355 x 355 x 355	45.0	N/A	N/A
HAO -3	455 x 455 x 455	95.0	N/A	N/A
HAO -4	455 x 455 x 605	125.0	N/A	N/A
HAO -5	605 x 605 x 605	224.0	N/A	N/A
HAO -6	605 x 455 x 910	252.0	N/A	N/A
HAO -7	605 x 605 x 910	336.0	N/A	N/A

Advance Quality Traders

HQ & Registered Office
Advance Quality Traders
D 14/31 Tedhi Neem Dashaswmedh Varanasi-221001
Tel: +91 0542-2419455 Fax: +91 0542-2419455
Mobial No. 08707286293 , 09451887577
Mail: info@advancescience.in
advancequalitytraders@gmail.com
Web: www.advancescience.in

Advance Quality Traders

R&D FACTORY WORK OFFICE ADDRESS:
Aarazi No.MI-271/5 Lane No. 10, Mahamanapuri
Colony Behind Karaudi, ITI college Karaudi Varanasi 221005
GST NO.:09BFKPS2845P1Z1
PHONE: 9451887577, 07355152295
EMAIL: info@advancescience.in
E-mail: info@advancescience.in
Web: www.advancescience.in